

TECHNICAL INFORMATION

Model/type	JKH-109Complex Compression Boot	Weight	1000g
Power supply	4500mAh lithium battery	Main unit size	165x83x57mm
Sleeve size	820x290mm	Degree of protection against electric shock	Type BF applied part
Rated voltage	3.7V	Type of protection against electric shock	Internally powered equipment (Not applicable)
Working current	600mA-2500mA	Grade of waterproof	IP22
Air pressure	0-120mmHg	Product life	1 year
Modes	2 modes	Software version	A0
Note: Not intended to be sterilized.			
Not for use in an OXYGEN RICH ENVIRONMENT			